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# Askina® Calgitrol® made easy

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Reference	Title	Type	Purpose	Outcome
<i>In vitro studies</i>				
Thomas S, McCubbin P. <i>J Wound Care</i> 2003; 12(8): 305-8.	Antimicrobial properties of ten silver-containing dressings	<i>In vitro</i>	To test 10 silver-containing dressings for the ability to produce sufficiently high concentrations of silver ions for antimicrobial efficacy	Calgitrol® Ag was shown to have sustained antimicrobial activity over two or more days Calgitrol® Ag was active within the first two hours of application, had no risk of transmission of microorganisms, and had best overall performance score
Aramwit P, et al. <i>Int J Molec Sci</i> 2010; 11: 2864-74	Evaluation of the antimicrobial effectiveness and moisture binding properties of wound dressings	<i>In vitro</i>	To investigate the antimicrobial efficacy and moisture penetration of five commercially available wound dressings	Calgitrol® Ag and ACTICOAT™ had the largest zone of inhibition, possibly due to the highest concentrations of silver All dressings, except Bactigras*, showed bactericidal activity, achieving a four log reduction in <i>E.coli</i> Calgitrol® Ag maintained the best environment within the wound with regards to moisture
<i>Clinical studies</i>				
Ricci E, et al. <i>Acta Vulnologica</i> 2007; 5(3): 105-11	Askina® Calgitrol® Ag: clinical use of an advanced ionic silver dressing	Open ended, non-controlled, non-randomised cohort study (n=37)	To determine the clinical efficacy of Calgitrol® Ag on patients with infected or critically colonised wounds; patients who required systemic antibiotics were not included in the study	Clinical signs of infection were alleviated in 34/37 (>90%) cases within the two week test period. No allergic reactions to the dressing occurred
Costa I, Linhares V. <i>Proceedings EWMA Conference 2008 (Lisbon)</i>	Second-degree burn in a diabetic – application of a silver alginate dressing	Case study with photo control	To evaluate the use of Calgitrol® Ag in a diabetic patient with a second degree burn on the arm	Complete healing in two weeks with a good aesthetic and functional result
Goeman C, et al. <i>Poster presentation. EWMA Conference 2008 (Lisbon)</i>	Dressing changes every 3 days: fiction or reality?	Case study with photo control	To examine whether the use of a silver alginate foam dressing with a change every 3 days is cost-effective	The use of silver alginate foam dressing with a renewal every 3 days is cheaper than a daily wound care with silver sulfadiazine ointment combined with an absorbent dressing
Durante CM. <i>Proceedings WUWHS Conference 2008 (Toronto)</i>	Clinical effectiveness of a polyurethane foam covered with a layer of calcium alginate and ionic silver	Prospective, open label (n=20)	To evaluate the clinical efficacy of Calgitrol® Ag in controlling the microbial development and managing the exudate	The use of Calgitrol® Ag lead to a decrease in the bacterial count and improvement of the local wound condition There was an absence of local and systemic complications No loss of absorption power under compression was found
Trial C, et al. <i>J Wound Care</i> 2010; 19(1)	Assessment of the antimicrobial effectiveness of a new silver alginate wound dressing: a RCT	Prospective, open-label, randomised controlled trial (n=42)	To compare the efficacy and tolerability of Calgitrol® Ag against Algosteril™ (a silver-free alginate dressing)	Calgitrol® Ag dressing had superior antimicrobial effect to Algosteril™ The two dressings were similar in terms of reduction of local infection, local tolerance, acceptability and usefulness
Opasonon et al. <i>Int Wound J</i> 2010; 7(6): 467-71	Clinical effectiveness of silver alginate dressing in outpatient management of partial-thickness burns	Prospective, descriptive study (n=65)	To compare the efficacy of the Calgitrol® Ag dressing and 1% silver sulfadiazine in patients with partial-thickness burn wounds	Average pain scores were lower ( $p<0.02$ ) in the Calgitrol® Ag group compared to the 1% silver sulfadiazine group The Calgitrol® group also had fewer dressing changes ( $p<0.02$ ), decreased nursing time ( $p<0.02$ ) and shorter healing time (7 vs 14 days in controls)
Chymrev IV. <i>Proceedings EWMA Conference 2011 (Brussels)</i>	Use of silver containing wound dressings after late necrectomy in the patients with deep burns	Cohort, prospective control (antiseptic ointment and polyurethane film) n=26 vs Calgitrol® Ag n=22	To compare the use of bandages plus antiseptic ointment and polyurethane film (control) with Calgitrol® Ag in patients with deep burns undergoing late necrectomy	By day 18 bacterial load was lower in Calgitrol® Ag treated patients than in the control group Wound preparation for split skin grafting was significantly shorter for patients treated with Calgitrol® Ag ( $20.6\pm1.9$ days vs $28.8\pm2.3$ days; $p<0.05$ )

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### To cite this publication

Opasanon S, Magnette A, Meuleneire F, Harding K. Askina® Calgitrol® Made Easy. *Wounds International* 2012; 3(1). Available from <http://www.woundsinternational.com>

## References

1. Leaper DJ. Silver dressings: their role in wound management. *Int Wound J* 2006; 3(4): 282-94.
2. Cutting K, White R, Edmonds M. The safety and efficacy of dressing with silver - addressing clinical concerns. *Int Wound J* 2007; 4(2): 177-84.
3. Thomas S, McCubbin P. Antimicrobial properties of ten silver-containing dressings. *J Wound Care* 2003; 12(8): 305-8.
4. BBraun 2011. Askina Calgitrol® Ag/Askina Calgitrol® THIN Ionic silver alginate dressings. BBraun Hospicare Ltd. Available at: <http://tinyurl.com/72rvwl>
5. Aramwit P, Muangman P, Namviriyachote N, Srichana T. In vitro evaluation of the antimicrobial effectiveness and moisture binding properties of wound dressings. *Int J Molec Sci* 2010; 11: 2864-74.
6. Goetz A. Water sanitation with silver. *J Am Water Works Assoc* 1943; 35: 579.
7. Rochat C, Uzdins K. Katadyn (silver preparation) clinical application. *Schweiz med Wochschr* 1947; 77: 1100-4.
8. Lansdown ABG. Silver 2: toxicity in mammals and how its products aid wound repair. *J Wound Care* 2002; 11(5): 173-77.
9. Jung WK, Koo H, Kim KW, et al. Antibacterial activity and mechanism of action of the silver ion in *Staphylococcus aureus* and *Escherichia coli*. *Appl Environ Microbiol* 2008; 74: 2171-78.
10. BBraun B 2011. Study of silver migration from wound dressing materials; Report HOSP217. BBraun Hospicare Ltd. Available from: <http://tinyurl.com/75gnnuc>
11. Edwards-Jones V. Antimicrobial and barrier effects of silver against methicillin-resistant *Staphylococcus aureus*. *J Wound Care* 2006; 15: 285-90.
12. Warriner R, Burrell R. Infection and the chronic wound. A focus on silver. *Adv Skin Wound Care* 2005; 18 Supp 1: 1-12.
13. BBraun 2011. Study of silver migration from wound dressing materials; Report HOSP216. BBraun Hospicare Ltd. Available from: <http://tinyurl.com/75gnnuc>
14. Roberts C, Ivins N, Widgerow A. ACTICOAT™ and ALLEVYN™ Ag Made Easy. *Wounds International* 2011; 2(2); Available from: <http://woundsinternational.com>
15. BBraun 2011. Study of silver migration from wound dressing materials; Report HOSP201. BBraun Hospicare Ltd. Available from: <http://tinyurl.com/7cqbh7>
16. BBraun 2011. Evaluation of leachable silver from a wound dressing using the swine model. Report HOSP213. BBraun Hospicare Ltd. Available from: <http://tinyurl.com/75gnnuc>
17. Trop M. Silver-coated dressing Acticoat caused raised liver enzymes and agryia-like symptoms in burn patient. *J Trauma* 2006; 61(1): 239-40.
18. Senior JR. Monitoring for hepatotoxicity: What is the predictive value of liver "function" tests? *Clin Pharmacol Ther* 2009; 85: 331-34.
19. BBraun 2011. Comparative evaluation of erythema and edema skin reaction and skin discoloration. Reports HOSP240. BBraun Hospicare Ltd. Available from: <http://tinyurl.com/75gnnuc>
20. Bradbury S, Fletcher J. Prontosan Made Easy. *Wounds International* 2011; 2(2). Available from <http://www.woundsinternational.com>
21. Opasanon S, Muangman P, Namviriyachote N. Clinical effectiveness of silver alginate dressing in outpatient management of partial-thickness burns. *Int Wound J* 2010; 7(6): 467-71.
22. Chymrev IV. Use of silver containing wound dressings after late necrectomy in the patients with deep burns. Proceedings EWMA Conference 2011, Brussels.
23. WUWHS. Wound infection in clinical practice. An international consensus. MEP Ltd: London, 2008.
24. Trial C, Darbas H, Lavigne JP. Assessment of the antimicrobial effectiveness of a new silver alginate wound dressing: a RCT. *J Wound Care* 2010; 19(1): 20-6.
25. Ricci E, Pittarello M, Cassino R, et al. Askina Calgitrol® Ag: clinical use of an advanced ionic silver dressing. *Acta Vulnologica* 2007; 5(3): 105-11.
26. Goeman C, Meuleneire F, Boucque H. Dressing changes every 3 days. Fiction or reality? Poster presentation. EWMA Conference, Lisbon 2008.